

Mallets

Martelos

Marteaux

Ø 30 mm, 300 g



1 plain face, 1 convex face
1 face plana e 1 face convexa
une face plane, l'autre convexe

COTTLE
PAN-41.19
19 cm - 7½"

The round face I used for the
Cartilage crusher PAN-34.70

Aface convexa se usa para prensa
De substância cartilaginosa PAN-34.70

La face convexe s' emploie pour presser la
Matière cartilagineuse PAN-34.70



QUISLING
PAN-42.17
17 cm - 6¾"

Intranasal mallet
Martelo intranasal
Marteau intranasal