




Wire forceps

Drahtzangen
Pinzas para hilos
Pincas à fil


Flat nose pliers, Wire twisting instruments

Flachzangen, Drahtspanninstrumente
Alicates planos, Tensores cierra-hilos
Pincas plates, Tendeurs serre-fils




PU. 80.01 PU. 80.02

17 cm




PU. 82.15
15 cm



PU. 83.18 PU. 84.18 PU. 85.18

18 cm



PU. 86.18



PU. 88.18 18 cm



PU. 89.17 17 cm



LOUTE
PU. 90.21 21 cm



Wire cutting pliers, PROLONG TC®
 Drahtkneifzangen, **PROLONG TC®**
 Alicates cortaambres, **PROLONG TC®**
 Pincas coupe-fil, **PROLONG TC®**

Pliers with tungsten carbide cutting edges

PU. 100.47
 For cutting off
 KIRCHNER wires,
 TIBIA head screw,
 bolts, pins upto
 ø 6.0 mm

PU. 100.56
PROLONG TC®
 For hard wires upto
 ø 2.5 mm



For hard wires upto
 ø 2.5 mm



PU. 98.22 22 cm

For hard wires upto
 ø 1.6 mm

For soft wires upto
 ø 2.0 mm



PU. 96.17 17 cm

MEGA CUT
PU. 100.47 47 cm
PU. 100.56 56 cm