

Instruments with Tungsten Carbide "TC"

## Professional "TC"



Non-glaring matt finish  
en mate deslumbrante  
exécution amte anti-éblouissante  
satinatura antiriflessi

Gold coloured  
for easy identification

Mangos dorados para  
mejor indentificación

bouts de ressorts dorés pour  
meilleure identification

parte inferiore dorata come  
identificazione speciale

Jaws fitted with Tungsten Carbide alloy, always guarantee an impeccable edge for sharpening as well as allowing for flexible and regular operations.

Bocas armadas con dos placas de cromo-tungsteno, una aleación dura comparable al diamante. Las estrías finas permiten coger los tejidos de forma segura y atraumática.

Plaques dures d'un alliage de chrome-tungstène avec mors dentelés fins pour une prise atraumatique et sûre des tissus sensibles.

Morsi con massiva placche in tungsteno, lega diamantata di cromo e wolframio, con denti molto rifiniti per una presa atraumatica sui tessuti delicati

Vacuum hardened surface for constant stability.

Temple de vacío en todas partes para flexibilidad y estabilidad constante

durci sous vacuum sur toute la surface pour élasticité et stabilité constante

tempra vacuum totale per una elasticità e stabilità permanente

Joined together to guarantee flex permanent flexibility

Resorte hecho con todo esmero garantizando una tensión suave y estable

ressort soigneusement travaillé garantissant une élasticité souple et permanente

molla finemente lavorata per una morbida elasticità permanente

## Atraumatic Forceps "TC"

Pinzas atraumáticas con placas de tungsteno  
Pincas atraumatiques avec plaquettes de tungstène  
Pinzette atraumatiche con morsi in tungsteno

### Instruments with Tungsten Carbide "TC"

Tissue holding and  
Suture Forceps


Pinzas para coger tejidos  
y para sutura

Pince pour saisir les tissus  
et à suture

Pinzette per tessuti  
e sutura



Adson  
**PF-15.62**  
12 cm-4 $\frac{3}{4}$ "  
"TC"







Adson  
**PF-15.65**  
15 cm-6"  
"TC"




Adson  
**PF-16.62**  
12 cm-4 $\frac{3}{4}$ "  
1 x 2 teeth  
"TC"

**Instruments with Tungsten Carbide "TC"**

**PF-01.64**

**PF-01.68**

**PF-01.70**

**PF-01.75**

**PF-39.65**

**PF-39.68**

**PF-39.70**

**PF-39.73**

**PF-39.75**
**Standard**
**PF-01.64** 14.5 cm-5½"

**PF-01.68** 18 cm-7"

**PF-01.70** 20 cm-8"

**PF-01.75** 25 cm-10"

**"TC"**
**Potts-Smith**
**PF-39.65** 15 cm-6"

**PF-39.68** 18 cm-7"

**PF-39.70** 20 cm-8"

**PF-39.73** 23 cm-9"

**PF-39.75** 25 cm-10"

**"TC"**

## Atraumatic Forceps "TC"

Pinzas atraumáticas con placas de tungsteno  
Pinces atraumatiques avec plaquettes de tungstène  
Pinzette atraumatiche con morsi in tungsteno

### Instruments with Tungsten Carbide "TC"



Dissector End  
Con cabo disector  
bout dissecteur  
terminale con  
reaschiatoio

Cushing  
**PF-52.67**  
17.5 cm-6¾"  
straight•recta•droie•retta  
"TC"



Cushing  
**PF-52.68**  
18 cm-7"  
angled•angulada•coudé•angolata  
"TC"



Dissector End  
Con cabo disector  
bout dissecteur  
terminale con  
reaschiatoio

Cushing  
**PF-55.67**  
17.5cm/ 7"  
"TC"